



Glorifico™



1. Make sure the glorifico is fully charged before use (green light indicator).

Make sure the areas you want to treat are dry. Do not put on any lotion before use.

2. After you turn it on it will start heating up slowly (no need to wait), place the probe (metal head) on your skin, the red light will then automatically activate.

3. Apply light pressure while slowly moving the glorificio in small circular motions.

4. apply your favorite skin care and treat each area for 3 more minutes with the silicon cap: face, forehead, eyes and neck

5. It is recommended to apply a small amount of moisturizing cream on the treated skin areas after each session and to massage it gently in circular motions. Any cream you'll use after treatment will be more effective.

6. Clean the glorifico after every use with a dry cloth or paper towel. Make sure not to allow any moisture to penetrate the seams of the unit. Don't get it wet.

the glorifico comes with 3 years warranty